

THERAFLEX-ECP

Extracorporeal Photochemotherapy

UVA illumination EVA bag - Ref. XUV8501Q 

Specifications

Nominal volume	Volume	Connections
3 L	300 mL	1 spike 1 injection site



8-MOP solution - Ref. CSL1800XA



Specifications

Volume	Composition
5 mL	Chemical composition for 100 mL of solution:
	Methoxalen.....2 mg Ethanol (96%).....3230 mg Glacial acetic acid.....230 mg Sodium acetate trihydrate.....860 mg Sodium chloride.....40 mg Water for injections.....100 mL

Produit Thérapeutique Annexe (PTA)
Marketing Authorisation in France

AMM PTA N°: TA 07.10.109

Macogenic G2 UVA illumination device - Ref. 9MG2000 

Specifications

Dimensions	Included items and options	Weight
87 cm (width) x 64 cm (depth) x 37 cm (height)	Printers, scanner, LCD colour touch screen	65 kg

