



DOCON CE1275

State-of-the-art blood collection mixer

**DESIGNED
TO EASE
YOUR WORK**



BENEFITS

- Large bag tray, smart tube detector, large display, remote control with adjustable height.
- Long life batteries allowing up to 1 week or 70 donations without recharge.
- High weighing precision (< 1%) and efficient mixing.
- High level of traceability through a large variety of data exchange features.
- Fully customizable collection programs.
- Integrated tube sealer (MacoSeal) as an option.



SEALING FUNCTION/MACOSEAL

Maintains the integrity of the donation. The Docon is equipped with the MACOSEAL quality hand sealer. Sealing operations can be recorded in the Docon with the donation data.



Docon[®], *Donation Controller*, is a **FULLY AUTOMATIC** blood collection mixer that guarantees **OPTIMAL** blood collection and maintains the requirements for **HIGH QUALITY** blood components.

Large bag tray to receive all types of blood bags.

Very large display (9 cm x 12 cm) :

information on the current status of the donation (duration of the donation, date, time and battery status). Measurement of blood flow and a flow graphic can be displayed upon demand.

Sensoclamp

Equipped with a smart tube recognition to guarantee the highest security during collection.

Automatic clamping upon insertion of the tube.



Robust transport case.

Docon Mobile II

Reduces the need to bend down. Shows the current status of the donation.

Provided with a keyboard : manual input of information such as temperature, haemoglobin etc.

Can be operated from a flexible gooseneck holder.



- Very quick start-up
- Fully customized and adapted to each specific need. Ex: type of bags, different sounds and visual alarms, bar code controls.

Go to www.macopharma.com to get a complete overview of our products.

DOCON is CE certified, Class IIa MDD, 93/42/EEC.

MACO PHARMA
Rue Lorthiois
F-59420 Mouvaux (France)
www.macopharma.com

Distributed by:

