




PLATELET CONCENTRATE FILTRATION

Ref. TRV8006XU



SPECIFICATIONS

Technical specifications

Product code	Bag(s)	Volume	Filter	Packaging	
TRV8006XU	2 PVC-DEHP 1 PVC-BTHC (Oxy-PL™)	1800 mL	THROMBOFLEX TXP	1 device / peelable sachet	15 devices / box

Included items

- SCD compatible tubings
- Tubing dedicated for the Additive Solution Bag
- Screw on / off sampling site

Sterilization

Steam

Instructions for use

NOTUNIPLT

Use

Leucocyte removal from Platelet Concentrate

Label

Universal with pictograms (according to ISO 3826-2) – ISBT128

Shelf Life

1 year



FINAL PRODUCT

- Leucoreduced Platelet Concentrate
350 mL +/- 10 %

TRV8006XU_A_A 24/02/2014