




COLLECTION BAGS 0459

Ref. MRT6280LS



SPECIFICATIONS

Technical specifications

| Product code | Bag(s) | Collection volume (mL) | Packaging | |
|--------------|------------|------------------------|----------------------------|------------------|
| MRT6280LS | 3 PVC-DEHP | 450 mL +/- 10 % | 1 device / peelable sachet | 15 devices / box |

Included items

- Needle (16G)
- Bactivam, Secuvam and Vacuvam
- Break-Away Cannula in the Collection Line

Solution

63 mL CPD - 100 mL SAGM

Sterilisation

Steam

Label


Spanish – ISBT128

Shelf Life

2 years

Instructions for use

NOTEURO04



FINAL PRODUCTS

- Plasma
- Buffy-Coat
- Red Cell Concentrate